SHIELDED BY SCIENCE : PIPELINE SHIELD STUDY ACROSS 204 PATIENTS.

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Hemorrhagic Stroke
Original research

Periprocedural to 1-year safety and efficacy outcomes with the Pipeline Embolization Device with Shield technology for intracranial aneurysms: a prospective, post-market, multi-center study


Hal Rice ¹, Mario Martínez Galdámez ², Markus Holtmannspötter ³, Laurent Spelle ⁴, Konstantinos Lagios ⁵,  Maria Ruggiero ⁶, Pedro Vega ⁷, Hemant Sonwalkar ⁸, René Chapot ⁹, Saleh Lamin ¹⁰

Correspondence to Dr Hal Rice, Interventional Neuroradiology, Gold Coast University Hospital, Southport, QLD 4215, Australia;

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
PRIMARY AND SECONDARY ENDPOINTS

PRIMARY EFFICACY ENDPOINT	Complete aneurysm occlusion defined as Raymond-Roy grade 1 without significant parent artery stenosis (<50%) or retreatment of the target aneurysm through 1 year.
PRIMARY SAFETY ENDPOINT	Major stroke in the territory supplied by the treated artery or <u>neurologic death</u> at 1-year post-procedure
Stroke is defined as a focal neurological deficit of presumed vascular origin persisting for >24 hours from symptom onset with a neuroimaging study or other quantitative study	
SECONDARY ENDPOINT #1	Major stroke in the territory supplied by the treated artery or neurological death at 30 days post-procedure due to procedural complications
SECONDARY ENDPOINT #2	Delayed intracerebral hemorrhage (ICH) >30 days post-procedure
SECONDARY ENDPOINT #3	Device deployment success rate at the target site

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Through 1 year post-procedure

- ✓ **Primary effectiveness endpoint: 71.7%**
 - **Complete Aneurysm Occlusion** in 75.0% of patients
 - **No significant stenosis (≥50%)** in 98.9% of patients
- ✓ **Primary Safety endpoint: 3.23%**
 - **Major stroke** experienced by 2.9% of patients
 - **Neurological death** occurred in 2 cases (1.0%)

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SECONDARY ENDPOINT #1 (SAFETY)

Major stroke in the territory supplied by the treated artery or neurological death at 30 days post-procedure due to procedural complications:
2.9% (6/204)

SECONDARY ENDPOINT #2 (SAFETY)

Delayed intracerebral hemorrhage (ICH) >30 days post-procedure:
0% (0/204)

SECONDARY ENDPOINT #3 (EFFICACY)

Device deployment success rate at the target site
98.0% (200/204)

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SHIELDED BY SCIENCE.

Pipeline™ Flex
Embolization Device
with Shield Technology™

SHIELD Study Results Across 204 Patients

REDUCES THROMBUS FORMATION.

3.2%

Primary Safety
Endpoint

PROMOTES HEALING.

77.2%

Complete Aneurysm
Occlusion at 12 Months

IMPROVES PERFORMANCE.

93.1%

Complete Wall Apposition
Post-Procedure

Hemorrhagic stroke

ORIGINAL RESEARCH

Periprocedural to 1-year safety and efficacy outcomes with the Pipeline Embolization Device with Shield technology for intracranial aneurysms: a prospective, post-market, multi-center study

Hal Rice,¹ Mario Martínez Galdámez,² Markus Holtmannspötter,³ Laurent Spelle,⁴ Konstantinos Lagios,⁵ Maria Ruggiero,⁶ Pedro Vega,⁷ Hemant Sonwalkar,⁸ René Chapot,⁹ Saleh Lamin¹⁰

ABSTRACT
Background The first and second generations of the Pipeline Embolization Device (PED) have been widely used for the treatment of intracranial aneurysms. Flow diversion addresses several of these limitations and is now broadly adopted as a minimally-invasive approach to conventional open microsurgical treatment of intracranial aneurysms.

J Neurointerv Surg: first published as 10.1136/neurintsurg-2017-013111

Click here to
download the
full article >