EMBOTRAP® II

DESIGNED TO COMPRESS LESS,  
TRAP MORE

Original Contribution

Primary Results of the Multicenter ARISE II Study  
(Analysis of Revascularization in Ischemic  
Stroke With EmboTrap)

Osama O. Zaidat, MD, MS\*; Hormozd Bozorgchami, MD; Marc Ribó, MD;  
Jeffrey L. Saver, MD; Heinrich P. Mattle, MD; René Chapot, MD; Ana Paula Narata, MD;  
Olivier Francois, MD; Ashutosh P. Jadhav, MD, PhD; Jonathan A. Grossberg, MD;  
Christian H. Riedel, MD; Alejandro Tomasello, MD; Wayne M. Clark, MD;  
Hannes Nordmeyer, MD; Eugene Lin, MD; Raul G. Nogueira, MD; Albert J. Yoo, MD, PhD;  
Tudor G. Jovin, MD; Adnan H. Siddiqui, MD, PhD; Thierry Bernard, MS; Mairsil Claffey;  
Tommy Andersson, MD, PhD\*

DESIGNED TO COMPRESS LESS,  
TRAP MORE

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ARISE II (Analysis of Revascularization in Ischemic Stroke With EmboTrap)는

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
objective	To evaluate EmboTrap's safety and efficacy compared with established stent retrievers, SWIFT and TREVO2 Trials
design	Single arm study Prospective clinical study Independent imaging core lab, CEC and DSMB
term	Oct 2015 - Feb 2017
hospitals	US 11, Europe 8
<b>Primary End Point</b>	<u>efficacy</u> : mTICI $\geq$ 2b within 3 passes of ET2 without the use of rescue therapy <u>Safety</u> : sICH within 24 hours Serious adverse device effect

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
- mTICI 3: 52%
- mTICI 2c-3: 76%
- mTICI 2b-3: 92.5%

#### Final TICI 2b-3\*


\*\*ARISE II study reports the end of procedure revascularization results which includes patients who had been treated with other therapies post initial use of EMBOTRAP II Device. Direct comparisons cannot be made.

<sup>^</sup>Outcomes generated from a Randomized Controlled Trial.

Berkhemer OA, et al. A randomized trial of intraarterial treatment for acute ischemic stroke. N Engl J Med 2015. (MRCLEAN)  
 Campbell BC, et al. Endovascular therapy for ischemic stroke with perfusion-imaging selection. N Engl J Med 2015. (EXTEND IA)  
 Goyal M, et al. Randomized assessment of rapid endovascular treatment of ischemic stroke. N Engl J Med 2015. (ESCAPE)  
 Saver JL, et al. Stent-retriever thrombectomy after intravenous t-PA vs. t-PA alone in stroke. N Engl J Med 2015. (SWIFT PRIME)  
 Jovin TG, et al. Thrombectomy within 8 hours after symptom onset in ischemic stroke. N Engl J Med 2015. (REVASCAT)  
 Zaidat OO, et al. Primary Results of the Multicenter ARISE II Study (Analysis of Revascularization in Ischemic Stroke With EmboTrap). Stroke 2018. (ARISE II)

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
#### First Pass mTICI

- One pass mTICI 3 – 30%
- One Pass mTICI 2c-3 – 40.1%
- One Pass mTICI 2b-3 – **51.5%**

#### Clinical Outcome

	ARISE2	SEER (1)	HERMES	REVASCAT*	SWIFT PRIME*	MR CLEAN*	DAWN*	EXTEND IA*	ESCAPE* (4)	THRACE*	ASTER-SR*	ASTER-ASP*
mRS 0-2 (5)	67.3%	56%	46%	44%	60%	33%	49%	71%	52%	53%	50%	45%

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

EMBOTRAP® II Revascularization Device does **not only** rely on radial force to extract clot, **COMPRESS LESS TRAP MORE**

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**Innovative, Dual Layer Design to ENGAGE & GRIP differently**

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Catalog #		Recommended Vessel Diameter	Microcatheter compatibility	Working Length	Total Length	Device Length	Tip Length
ET009521	5 x 21	1.5 - 5.0mm	0.021" ID	21mm	33mm	194cm	4mm
ET009533	5 x 33			33mm	44mm	195cm	4mm