



The PERFECT study aimed to characterize the performance of EMBOVAC™ in the treatment of AIS in a real-world clinical setting.

PERFECT 연구는 실제 임상 환경에서 EMBOVAC™의 AIS 치료 효과를 파악하기 위해 시행되었습니다.

## METHODOLOGY

Prospective, multi-center, single arm, post-market clinical follow-up study

October 2020 — July 2022



100

AIS PATIENTS WITH LVOS



11

EUROPEAN CENTERS

### TECHNIQUE



**Direct aspiration with EMBOVAC™ was the first-line therapy.** A minimum of 3 passes with EMBOVAC™ were required before switching strategy.

### DATA ASSESSMENTS



Imaging and procedure angiography were assessed by an independent core lab



90-day mRS assessed by an independent evaluator



Device related adverse events and symptomatic ICH assessed by an independent clinical events committee



Retrieved clot was collected per pass from 83 subjects, and clot was analyzed by an independent central clot lab

### OUTCOMES



**PRIMARY OUTCOME**  
Successful reperfusion at the end of procedure  
Final mTICI  $\geq 2b$



**SECONDARY OUTCOME**  
• First pass mTICI  $\geq 2c$   
• mRS 0-2 at 90 days  
• 90-day all-cause mortality

AIS, acute ischemic stroke; LVO, large vessel occlusion; mRS, modified Rankin Scale; mTICI, modified Thrombolysis in Cerebrovascular Infarction; sICH, symptomatic intracranial hemorrhage

Lobotesis, K. et al, Direct Aspiration with EMBOVAC — First Clinical Experience and clot Composition (Results of the "PERFECT" Study. In ESMINT (2023).

# RESULTS

## PATIENT DEMOGRAPHICS AND BASELINE CHARACTERISTICS

|                               |                |
|-------------------------------|----------------|
| Age                           | 70.4±13.95     |
| Sex, F                        | 59.0% (59/100) |
| Pre-stroke mRS 0-1            | 98.0% (98/100) |
| Baseline NIHSS                | 14.9±6.39      |
| Use of IV-tPA                 | 51.0% (51/100) |
| Baseline ASPECT*              |                |
| 0-5                           | 15.0% (15/100) |
| 6-7                           | 25.0% (25/100) |
| 8-10                          | 60.0% (60/100) |
| Witnessed stroke              | 72.0% (72/100) |
| Wake-up stroke                | 9.0% (9/100)   |
| Unwitnessed non-wakeup stroke | 19.0% (19/100) |
| Occlusion Location*           |                |
| ICA & Carotid T               | 18.0% (18/100) |
| MCA, M1                       | 71.0% (71/100) |
| MCA, M2                       | 6.0% (6/100)   |
| Other                         | 4.0% (4/100)   |
| Cannot determine              | 1.0% (1/100)   |

\*core lab assessment

## PROCEDURAL OUTCOMES

**27.8 min**

Procedure time, arterial  
puncture to mTICI ≥ 2b

**2.4**

Average number  
of total passes

**2.0**

Average number of  
EMBOVAC™ passes

ASPECT, Alberta Stroke Program Early CT Score; ICA, internal carotid artery; IV-tPA, intravenous tissue-type plasminogen activator; MCA, middle cerebral artery; mTICI, modified Thrombolysis in Cerebrovascular Infarction; mRS, modified Rankin Scale; NIHSS, National Institute of Health Stroke Scale


Lobotesis, K. et al, Direct Aspiration with EMBOVAC — First Clinical Experience and clot Composition (Results of the "PERFECT" Study. In ESMINT (2023).

## RESULTS


The PERFECT study demonstrated direct aspiration with EMBOVAC™ as first line therapy resulted in high rates of reperfusion, good clinical outcomes, and a strong safety profile.

PERFECT 연구는 EMBOVAC™을 사용한 direct aspiration이 first line 치료로 시행될 경우 높은 혈류재개통, 긍정적인 임상 결과 및 안전성을 입증하였습니다.

### ANGIOGRAPHIC OUTCOMES



### CLINICAL OUTCOMES



### SAFETY OUTCOMES



<sup>a</sup> Estimate based on Kaplan-Meier analysis

mRS, modified Rankin Scale; mTICI, modified Thrombolysis in Cerebrovascular Infarction; SAE, serious adverse event; sICH, symptomatic intracranial hemorrhage

Lobotesis, K. et al., Direct Aspiration with EMBOVAC — First Clinical Experience and clot Composition (Results of the "PERFECT" Study). In ESMINT (2023).