Medtronic

Power to React

React & Solitaire - Clinical Evidence

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INSPIRE-S


Innovative Neurovascular Product Surveillance **RE**gistry-Acute Ischemic **Stroke**

A neurovascular registry to study innovative AIS approved devices

Study design

- Prospective, global registry¹
- >800 AIS patients enrolled to date, across 33 sites worldwide² (29 Western European countries)²
- Core-lab adjudicated
- Clinical events committee (CEC) adjudicated
- 3 cohorts based on first pass technique:

Purpose

To present the results for **397** patients from the combination therapy cohort.

Primary objective

Good functional outcomes as measured by a modified Rankin Scale (mRS) score of 0-2, or return to pre-stroke mRS at 90 days post-intervention.

Additional outcomes

Final successful revascularization, first pass revascularization, and safety outcomes.

Procedural characteristics

Combination Therapy % (n/Total N)


Thrombectomy Method at 1st Pass	
Solitaire™ (Platinum and X)	100% (397/397)
React™	76.1% (302/397)
React™ 68	37.8% (150/397)
React™ 71	38.3% (152/397)
Other Aspiration Catheter*	23.9% (95/397)
Balloon Guide Catheter used at 1st Pass	39.5% (157/397)

*Enrollment guide was updated in March 2022 to include use of other market-approved aspiration catheters in conjunction with Solitaire™ as the primary treatment device.


Preliminary Results

First Pass Revascularization


eTICI ≥2b50


eTICI ≥2b67

eTICI ≥2c

eTICI = 3

Final Successful Revascularization

Combination Therapy %

eTICI ≥ 2b50	93.4%
eTICI ≥ 2b67	88.8%
eTICI ≥ 2c	75.6%
eTICI = 3	52.7%

397 Patients treated with Combination Therapy

100% Solitaire™ (Platinum and X) used at 1st pass

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Conclusion

This analysis of INSPIRE-S supports the real-world safety and effectiveness of the use of the Solitaire™ Revascularization Device in Combination Therapy, primarily with the React™ Catheter, in patients with AIS.

INSPIRE-S

A neurovascular registry to study innovative AIS approved devices

References

¹INSPIRE-S Analysis Plan V3.0.

²Medtronic data on file for INSPIRE-S Registry

Presented by Marc Ribo at the European Stroke Conference Basel, 15-17 May 2024

Presented by Marc Ribo at ISC Phoenix, 7-9 February 2024

ISC <https://professional.heart.org/en/meetings/international-stroke-conference>

See the device manual for detailed information regarding the instructions for use, indications, contraindications, warnings, precautions, and potential adverse events. For further information, contact your local Medtronic representative and/or consult the Medtronic website at

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